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# IgA

## Immunoturbidimetric method for the determination of immunoglobulin A

### SUMMARY

The quantitative determination of IgA is needed for typify immunodeficiencies and myelomas, as well as for the follow-up of acute and chronic infections.

### PRINCIPLE

The immunoglobulin A reacts with specific antibody generating insoluble immune complexes. The turbidity produced by these immune complexes is proportional to the IgA concentration in the sample and can be measured spectrophotometrically.

### PROVIDED REAGENTS

**A. Reagent A:** buffered saline solution, pH 7.5.

**B. Reagent B:** antibody monospecific anti-IgA.

### NON-PROVIDED REAGENTS

- Saline solution
- Wiener lab.'s **Calibrador Proteínas Nivel Alto Turbitest AA**.

### INSTRUCTIONS FOR USE

**Provided Reagents:** ready to use.

### WARNINGS

Reagents are for "in vitro" diagnostic use.

All samples from patients should be handled as capable of transmitting infection.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

### STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable at 2-10°C until the expiration date stated on the box. Do not freeze.

### SAMPLE

Serum or heparinized plasma.

**a) Collection:** obtain in the usual way.

**b) Additives:** if plasma is used, it is recommended to use heparin as anticoagulant.

**c) Known interfering substances:**

- Do not use contaminated or hemolyzed sera.
- No interferences have been observed with bilirubin up to 20 mg/dl, triglycerides up to 25 g/l and hemoglobin up to 1 g/dl. See Young, D.S. in References for effect of drugs on the present method.

**d) Stability and storage instructions:** sample should be preferably fresh. In case the test cannot be performed on the day, the sample can be store for up to 1 week at 2-10°C. In case it must be processed later than that keep it at -20°C.

### REQUIRED MATERIAL (non-provided)

- Spectrophotometer.
- Spectrophotometric cuvettes.
- Micropipettes and pipettes for measuring the stated volumes
- Kahn or hemolysis tubes.
- Stopwatch.

### ASSAY CONDITIONS

- Wavelength: 340 nm
- Reaction temperature: room temperature (25°C). Temperature control is not critical, it can range between 22 and 30°C.
- Reaction time: 30 minutes

### PROCEDURE

#### CALIBRATION CURVE

In Kahn tubes dilute the Calibrador Proteínas Nivel Alto with saline solution 1:10, 1:20, 1:40, 1:80 and 1:160, using saline solution as the zero point.

<b>Diluted Calibrador Proteínas</b>	50 ul
-------------------------------------	-------

<b>Reagent A</b>	900 ul
------------------	--------

Homogenize and measure the absorbance of each dilution at 340 nm ( $OD_1$ ), setting the instrument to zero with distilled water. Then, add:

<b>Reagent B</b>	80 ul
------------------	-------

Mix and incubate 30 minutes at room temperature. measure the absorbance at 340 nm ( $OD_2$ ), setting the instrument to zero with distilled water.

Calculate the absorbance difference ( $\Delta A = OD_2 - OD_1$ ) for each Calibrador Proteínas dilution, including the zero point.

Draw on graph paper the  $\Delta A$  absorbance differences based on the Calibrador Proteínas concentration in mg/dl (g/l).

#### SAMPLES PROCEDURE

Dilute the samples 1:10 with saline solution.

<b>Diluted Sample</b>	50 ul
-----------------------	-------

<b>Reagent A</b>	900 ul
------------------	--------

Homogenize and measure the absorbance at 340 nm (OD<sub>1</sub>), setting the instrument to zero with distilled water. Then add:

<b>Reagent B</b>	80 ul
------------------	-------

Mix and incubate 30 minutes at room temperature. measure the absorbance at 340 nm (OD<sub>2</sub>), setting the instrument to zero with distilled water.

### CALCULATIONS

Calculate the absorbance difference ( $\Delta A = OD_2 - OD_1$ ) for each sample tested. Interpolate this  $\Delta A$  in the calibration curve to determine the concentration in mg/dl (g/l) corresponding to the sample under study. Samples with absorbance values higher than the absorbance measurement for **Calibrador Proteínas nivel alto** should be diluted 1:2 with saline solution and retested. Multiply the obtained result by 2.

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The Control should be processed in the same manner as samples.

### REFERENCE VALUES

Children and young people:

0-1 year: 0 - 83 mg/dl (0 - 0,83 g/l)

1-3 years: 20 - 100 mg/dl (0,20 - 1,00 g/l)

4-6 years: 27 - 195 mg/dl (0,27 - 1,95 g/l)

7-9 years: 34 - 305 mg/dl (0,34 - 3,05 g/l)

10-11 years: 53 - 204 mg/dl (0,53 - 2,04 g/l)

12-13 years: 58 - 359 mg/dl (0,58 - 3,59 g/l)

14-15 years: 47 - 249 mg/dl (0,47 - 2,49 g/l)

16-19 years 61 - 348 mg/dl (0,61 - 3,48 g/l)

Adults:

70 - 400 mg/dl (0.7 - 4.0 g/l)

Each laboratory should set its own reference values.

### SI SYSTEM UNITS CONVERSION

IgA (mg/dl) x 10 = IgA (mg/l)

### PROCEDURE LIMITATIONS

If the samples have no definite clinical characteristics, a protein electrophoresis should be performed to detect an unexpected antigen excess, as it happens in gammopathies. Turbidity and particles in the sample may interfere with the test. Therefore, the particles that may be the result of an incomplete coagulation or protein denaturalization must be removed by centrifugation before testing the sample.

### PERFORMANCE

**a) Reproducibility:** simultaneously processing 20 replicates of one sample, the following results were obtained:

Level	S.D.	C.V.
37 mg/dl	± 1.12 mg/dl	3.04%
192 mg/dl	± 4.92 mg/dl	2.56%
560 mg/dl	± 19.7 mg/dl	3.52%

**b) Dynamic range:** values can be obtained between the lowest and highest calibrator concentrations of the calibration curve (approximately 600 mg/dl).

**c) Detection limit:** the minimum detectable concentration change of IgA is 18 mg/dl.

### PARAMETERS FOR AUTOANALYZERS

Refer to the specific applications of each autoanalyzer. For calibration must be use Wiener lab's **Calibrador Proteínas nivel alto Turbitest AA** following the autoanalyzer requirements.

### WIENER LAB. PROVIDES

1 x 60 ml Reagent A

1 x 5 ml Reagent B

(Cat. 1513261)

1 x 60 ml Reagent A

1 x 5 ml Reagent B

(Cat. 1009344)

1 x 60 ml Reagent A

1 x 5 ml Reagent B

(Cat. 1009217)

1 x 60 ml Reagent A

1 x 5 ml Reagent B

(Cat. 1009667)

1 x 60 ml Reagent A

1 x 5 ml Reagent B

(Cat. 1009961)

2 x 36 ml Reagent A

2 x 3 ml Reagent B

(Cat. 1008132)

### REFERENCES

- Ichihara, K. et al - J. Clin. Lab. Anal. 10:110, 1996.

- Itoh, Y. et al - J. Clin. Lab. Anal. 11:39,1997.

- Maynard, Y. et al - Clin. Chem. 32/5:752, 1986.

- Dati, F. - Journal of IFCC VIII/1:29, 1996.

- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4<sup>th</sup> ed., 2001.

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Wiener Laboratorios S.A.I.C.